

CLEAN LINEN TROLLEY 2 doors - medium



Product Code	45933
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279459333

Description

CLEAN LINEN TROLLEYS 2 DOORS

Clean linen trolley - medium
Size (mm): 820x640xh 1280

Solid and sturdy trolleys for distribution of clean linen and patient hygiene.
Steel painted with antimicrobial epoxy powders to obtain maximum lightness.
Upper part has large support surface in technopolymer plastic material.
Integrated ergonomic push handle.

Supplied with 2 doors with 3 internal shelves. Key closures. 4 castors Ø 125 mm, anti-wire, pivoting, two of which with lever locking system and reinforced plastic support.

Made in Italy.