



Product Code	45950
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8683907116900

Description

BASIC CLEANING TROLLEY

Basic trolley suitable for floor cleaning in large areas.

Equipped with:

- 2x 25 l buckets, blue and red, which allow to separate cleaning solution from rinse water
- wringer
- 4 swivelling castors Ø 80 mm, with antishock system

Size: 77x42x85 cm

Weight: 8.1 kg

Packing size: 80x43x58 cm

MULTILINGUAL INSTRUCTIONS: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic