



Product Code	45954
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8683907116986

Description

SAFETY MULTIFUNCTION CLEANING TROLLEY

Service trolley with pre-impregnated system, particularly suitable in the Ho.Re.Ca and hospital world. Sturdy and easy to clean.

Equipped with:

- 2 doors to store items from right and left side, with 3 shelves inside
- 2x 60 l double canvas bag holder for waste sorting
- 120 l canvas bag cover
- 2x 22 l containers for pre-soaked cloths, in blue and red colours
- 4x 6 l buckets with lids, in blue, red, yellow and green colours
- 3 broom hooks
- mop handle holder and frame holder

Packing size: 94x60x62 cm

Size: 122x57x119 cm

Weight: 29.8 kg

4 wheels Ø 100mm

MULTILINGUAL INSTRUCTIONS: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic