

LAERDAL STIFNECK SELECT EXTRICATION COLLAR - pediatric



Product Code	34628
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2239719
CND	Y060312
EAN/UPC	7045432105599

Description

LAERDAL STIFNECK SELECT EXTRICATION COLLAR - pediatric

This product is intended to support cervical spine in a neutral position during transportation, in combination with other cervical and full body immobilization devices (CSIDs and FBSIDs).

The Laerdal Stifneck select extrication collar is four different collars in one.

Stabilize spinal patients correctly and easily with 4 different sizing options that all lock quickly and safely into place.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, PL, HR, RO, RU, CZ, SI, GR, EE.