

PREVIMEDICAL COMPRESSION BANDAGE 5 m x 8 cm



Product Code	34654
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2371616
CND	M0304010103
EAN/UPC	8033140281604
PARAF	977814235

Description

ELASTIC COMPRESSION BANDAGES - individually packed, with clip.

PREVIMEDICAL - medium light compression

- universal compression bandage
- protects and supports all types of medication
- to fix medications, for support and relief in case of contusions and haematomas with medium-light compression

Technical Specifications

- Composition:
 - Cotton 60%
 - Polyamide 33%
 - Polyurethane 7%
- Extensibility:
 - Longitudinal 160%
 - Transversal -