

ELASTIC COHESIVE BANDAGE 4 m x 6 cm



Product Code	34657
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2371622
CND	M0303010201
EAN/UPC	8033140280690
PARAF	977814250

Description

- Elastic cohesive bandage:
- To fix dressings
 - Ideal for supporting light joint fractures
 - This bandage does not require clips because perfectly adherent

Technical Specifications

- Composition:
 - 42% cotton
 - 29% viscose
 - 29% polyamide, latex
- Width: 6 cm
- Length: stretched length 4 m
- Extensibility: approx. 50%