

COHESIVE COLOUR BANDAGE 4 m x 6 cm - mixed colours



Product Code	34668
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2269546
CND	M0303010201
EAN/UPC	8033140281789
PARAF	977814324

Description

- Cohesive bandages through microfine latex dots:
- For fixing dressings
 - Ideal for supporting joint treating light fractures
 - This bandage does not require clips because it remains well-attached to itself

Technical Specifications

- Composition:
 - 42% cotton
 - 29% viscose
 - 28% polyamide
 - 1% elastan and latex
- Width: 6 cm
- Length: stretched length 4 m
- Extensibility: approx. 90%