

# COHESIVE COLOUR BANDAGE 4 m x 8 cm - mixed colours



Product Code	34673
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2269546
CND	M0303010201
EAN/UPC	8033140281796
PARAF	977814375

## Description

- Cohesive bandages through microfine latex dots:
- For fixing dressings
  - Ideal for supporting joint treating light fractures
  - This bandage does not require clips because it remains well-attached to itself

## Technical Specifications

- Composition:
  - 42% cotton
  - 29% viscose
  - 28% polyamide
  - 1% elastan and latex
- Width: 8 cm
- Length: stretched length 4 m
- Extensibility: approx. 90%