



| | |
|---------------|----------------|
| Product Code | 34676 |
| Unit of sale | box of 10 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2269546 |
| CND | M0303010201 |
| EAN/UPC | 8033140281710 |
| PARAF | 977814401 |

Description

- Cohesive bandages through microfine latex dots:
- For fixing dressings
 - Ideal for supporting joint treating light fractures
 - This bandage does not require clips because it remains well-attached to itself

Technical Specifications

- Composition:
 - 42% cotton
 - 29% viscose
 - 28% polyamide
 - 1% elastan and latex
- Width: 10 cm
- Length: stretched length 4 m
- Extensibility: approx. 90%