



|               |                |
|---------------|----------------|
| Product Code  | 34682          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2528765       |
| CND        | M03050299     |
| EAN/UPC    | 8023279346824 |
| PARAF      | 903163590     |

## Description

These supports immobilize your limbs in an efficient and practical way.

Each splint can be shaped to fit the patient's personal need.

It ensures a high degree of rigid support: thanks to inner part made of neoprene lined with 100% laminated, absorbent-proof nylon assembled to a rigid metal splint.

A velcro band allows a quick fixing.

Available in 5 sizes or in a bag set of 5 pieces.

Instructions in: GB, FR, IT, ES, PT.