



Product Code	34684
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2528765
CND	M03050299
EAN/UPC	8023279346848
PARAF	903163614

## Description

These supports immobilize your limbs in an efficient and practical way.

Each splint can be shaped to fit the patient's personal need.

It ensures a high degree of rigid support: thanks to inner part made of neoprene lined with 100% laminated, absorbent-proof nylon assembled to a rigid metal splint.

A velcro band allows a quick fixing.

Available in 5 sizes or in a bag set of 5 pieces.

Instructions in: GB, FR, IT, ES, PT.