

GIMASPLINT - set of 5 pieces with bag



Product Code	34685
Unit of sale	1 set
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2528428
CND	M03050299
EAN/UPC	8023279346855
PARAF	903153385
GMDN	12095

Description

These supports immobilize your limbs in an efficient and practical way.

Each splint can be shaped to fit the patient's personal need.

It ensures a high degree of rigid support: thanks to inner part made of neoprene lined with 100% laminated, absorbent-proof nylon assembled to a rigid metal splint.

A velcro band allows a quick fixing.

Available in 5 sizes or in a bag set of 5 pieces.

Instructions in: GB, FR, IT, ES, PT.