

# VACUUM SPLINT KIT



Product Code	34690
Unit of sale	1 set
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2598644
CND	M03050201
EAN/UPC	6974442059183

## Description

### VACUUM SPLINT KIT

The kit includes vacuum splints for arm, elbow and leg, nylon bag and pump.  
Made of high quality and durable TPU material, filled with harmless EPS particulate material.  
Suitable for X-ray and CT scan.  
The firmness of the splints can be adjusted by pumping air according to the situation.  
Quick and easy softness and hardness adjustment of the splint.

Packaging size: 50x28x22 cm

Carrying bag and 0.25 l pump included