


SET OF LATEX TONE LOOPS



Product Code	47040
Unit of sale	set of 4 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279470406
PARAF	977817194

Description

SET OF 4 TONE LOOPS

The closed-loop flat band fits comfortable around the upper and lower legs. It helps to exercise the muscles of hips and legs.

Set includes 4 below loops: each in a different colour-coded level.

Tone loop 18 mm x 58 cm x 1.4 mm (X-Light - yellow)

Tone loop 24 mm x 58 cm x 1.4 mm (Medium - red)

Tone loop 24 mm x 58 cm x 1.6 mm (Heavy - green)

Tone loop 30 mm x 58 cm x 1.6 mm (X-Heavy - blue)