

LATEX EXERCISE TUBE WITH TPR HANDLES 125 cm x 1.5 mm - X-light - yellow



Product Code	47050
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8023279470505
---------	---------------

PARAF	977817206
-------	-----------

GTIN	08023279470505
------	----------------

Description

Latex resistive exercise tube with two pre-attached durable foam covered handles for comfort and grip. 5 colour-coded resistance levels available.

- Strength: yellow (X-Light)
- Length: 125 cm

Tubing is an affordable and convenient training tool which provides effective upper and lower-body workout.

Elastic properties for resistive and strengthening exercise.