

LATEX EXERCISE TUBE WITH TPR HANDLES 125 cm x 2.5 mm - medium - red



Product Code	47052
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279470529
PARAF	977817220
GTIN	08023279470529

Description

Latex resistive exercise tube with two pre-attached durable foam covered handles for comfort and grip. 5 colour-coded resistance levels available.

- Strength: red (Medium)
- Length: 125 cm

Tubing is an affordable and convenient training tool which provides effective upper and lower-body workout.

Elastic properties for resistive and strengthening exercise.