

# LATEX EXERCISE TUBE WITH TPR HANDLES 125 cm x 3.0 mm - heavy - blue



Product Code	47053
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8023279470536
---------	---------------

PARAF	977817232
-------	-----------

GTIN	08023279470536
------	----------------

## Description

Latex resistive exercise tube with two pre-attached durable foam covered handles for comfort and grip. 5 colour-coded resistance levels available.

- Strength: blue (Heavy)
- Length: 125 cm

Tubing is an affordable and convenient training tool which provides effective upper and lower-body workout.

Elastic properties for resistive and strengthening exercise.