



Product Code	47060
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8023279470604
---------	---------------

PARAF	977817257
-------	-----------

GTIN	08023279470604
------	----------------

## Description

Latex resistive exercise tube

- Thickness: 1.5 mm
- Strength: yellow (X-Light)
- Length: 25 m

Tubing is an affordable and convenient training tool which provides effective upper and lower-body workout.

Elastic properties for resistive and strengthening exercise.  
5 colour-coded resistance levels available.