



Product Code	47063
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279470635
PARAF	977817283
GTIN	08023279470635

Description

Latex resistive exercise tube

- Thickness: 3 mm
- Strength: blue (Heavy)
- Length: 25 m

Tubing is an affordable and convenient training tool which provides effective upper and lower-body workout.

Elastic properties for resistive and strengthening exercise.
5 colour-coded resistance levels available.