



Product Code	47064
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279470642
PARAF	977817295

## Description

Latex resistive exercise tube

- Thickness: 3.5 mm
- Strength: black (X-Heavy)
- Length: 25 m

Tubing is an affordable and convenient training tool which provides effective upper and lower-body workout.

Elastic properties for resistive and strengthening exercise.  
5 colour-coded resistance levels available.