

# BURST RESISTANT BALL diam. 65 cm - green



Product Code	47103
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279471038
---------	---------------

PARAF	977817358
-------	-----------

## Description

Designed for dynamic workouts that incorporate moderate level of resistance.

Deflate slowly if punctured.

136 Kg (300 lb) burst resistance level.