

ELASTICATED BANDAGES 5 cm x 4.5 m - box of 20



Product Code	34720
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2427635
CND	M0303010101
EAN/UPC	8023279347203
PARAF	977814591

Description

Elasticated non adhesive bandages for containment, individually cellophane wrapped, in a box.
Non adhesive bandages are made of 65% cotton and 35% polyamide fibre.
Provided with clasp.
Made in Italy.