

PREVICREPE FIXATION BANDAGES 4 m x 7 cm



Product Code	34726
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2371596
CND	M0303010101
EAN/UPC	8033140280287

Description

PREVICREPE S FIXATION ELASTICATED BANDAGES - 75 g2

For fixing dressing, for support and light compression.
Its characteristic crepe like touch enables the bandage to remain well-attached to itself.
Well-tolerated by the skin, it ensures a high level of breathing.
Washable at 60°C.
Made of 99% cotton - 1% lycra.
Made in Europe.