



Product Code	34727
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2371596
CND	M0303010101
EAN/UPC	8033140280294

## Description

### PREVICREPE S FIXATION ELASTICATED BANDAGES - 75 g2

For fixing dressing, for support and light compression.  
 Its characteristic crepe like touch enables the bandage to remain well-attached to itself.  
 Well-tolerated by the skin, it ensures a high level of breathing.  
 Washable at 60°C.  
 Made of 99% cotton - 1% lycra.  
 Made in Europe.