

COHESIVE NON WOVEN ELASTIC BANDAGE 4.5 m x 7.5 cm - red



Product Code	34735
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2269553
CND	M0303010201
EAN/UPC	8033140281871

Description

PREVITNT NON WOVEN ELASTIC BANDAGES

For fixing dressings. Ideal for supporting joints, treating light fractures.
This bandage does not require clips because it remains well-attached to itself.
Particularly suitable for use on joints, for sports and veterinary use.

- Composition:
 - 85% polypropylene
 - 15% spandex and latex.
- Extensibility: approx 100%.