

COHESIVE NON WOVEN ELASTIC BANDAGE 4.5 m x 10 cm - blue



|               |                |
|---------------|----------------|
| Product Code  | 34738          |
| Unit of sale  | box of 10 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2269553       |
| CND        | M0303010201   |
| EAN/UPC    | 8033140283004 |

Description

PREVITNT NON WOVEN ELASTIC BANDAGES  
For fixing dressings. Ideal for supporting joints, treating light fractures.  
This bandage does not require clips because it remains well-attached to itself.  
Particularly suitable for use on joints, for sports and veterinary use.

- Composition:
  - 85% polypropylene
  - 15% spandex and latex.
- Extensibility: approx 100%.