

COHESIVE NON WOVEN ELASTIC BANDAGE 4.5 m x 10 cm - 5 colours



| | |
|---------------|----------------|
| Product Code | 34740 |
| Unit of sale | box of 10 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2269553 |
| CND | M0303010201 |
| EAN/UPC | 8033140281918 |

Description

PREVITNT NON WOVEN ELASTIC BANDAGES

For fixing dressings. Ideal for supporting joints, treating light fractures.
This bandage does not require clips because it remains well-attached to itself.
Particularly suitable for use on joints, for sports and veterinary use.

Mix 2 bandages x 5 colours (red, green, blue, yellow, skin).

- Composition:
 - 85% polypropylene
 - 15% spandex and latex.
- Extensibility: approx 100%.