

KINESIOLOGY TAPE PRO 5 m x 5 cm - skin



Product Code	34757
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2450789
CND	M03040299
EAN/UPC	8023279347579
PARAF	989690944
GMDN	10284

Description

KINESIOLOGY PRO TAPES - thickness 0.32 mm

- Kinesiology PRO tapes 5 m x 5 cm - skin

High quality Kinesiology tape conforms to all areas of the body. Provides firm support and compression. Ideal for treatment of muscular symptoms, stimulus to blood circulation and lymph and fascia diversion. Reduce pain improving athletic performance. Easy to tear and unwinds to the core.

Individually packed in a plastic bag and cardboard box.

- soft hypoallergenic cotton cloth, elasticity up to 160%
- 70% water proof (level 3) and sweat proof
- supersticky: over 10 days

Kinesiology PRO: supersoft and quick drying

Multilingual box and instructions (GB, FR, IT, ES)

Made of 97% cotton, 3% spandex