



Product Code	44760
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2831159
CND	V0899
EAN/UPC	8023279447606
GTIN	08023279447606
GMDN	35892

Description

WARD STRETCHER

Main features:

- Platform: all-welded mainframe made of semi-oval steel tube 50x25x2 mm.
- Undercarriage: made of steel tube Ø 30x1.5 mm and oval steel tube 40x20x1.5 mm
- Push/Pull handle: made of steel tube Ø 30x1.5 mm.
- Surface finish: epoxy gloss finish standard grey 7040
- Mattress platform section: two section mattress platform fitted with grids made of 5 mm Ø steel rods and square mesh 50x50 mm.
- Backrest: square cross section 25x25x1.5 mm, adjustable by means of a rack device (5 positions).
- Castors - four castors Ø 200 mm made of grey synthetic rubber two of which with independent brake.
- Bumpers: Four bumpers at the corners made of anti-acid and anti-shock grey resin.

Supplied disassembled.

External size: 202x71xh 90 cm

Wired net height: 78 cm

Load: 135 kg.

Manual: GB, FR, IT, ES.

FLAT PACK reduced shipping volume

Made in Italy.