


OXYGEN BOTTLE HOLDER - chrome plated steel



Product Code	44764
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2278116
CND	V9099
EAN/UPC	8023279447644
GTIN	08023279447644
GMDN	35892

Description

Oxygen bottle holder - chrome plated steel (for oxygen bottle up to 41 cm, Ø 15,5 cm).

Accessoire for ward stretcher (code 44760)