

COTTON GAUZE SWABS 20x20 cm - bag of 25 pcs



34767

34766

34765

Product Code	34767
Unit of sale	100 bags
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	39269
CND	M0201010101
EAN/UPC	8000246004307
PARAF	931988063

Description

STERILE COTTON HYDROPHILIC GAUZE

- Size: 20x20 cm
- Threads: 12/8
- Minimum order: carton of 100 bags of 25 pcs

Folded hydrophilic gauze 100% pure cotton, suitable for sterile medical dressings, to buffer and protect wounds.
It allows skin perspiration.
Individually packed and latex free.

Instructions and box in GB, IT.
Made in Italy.