

COTTON GAUZE SWABS 36x40 cm - box of 12



Product Code	34769
Unit of sale	48 boxes
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	39903
CND	M0201010101
EAN/UPC	8000246004352
PARAF	931972069

Description

STERILE COTTON HYDROPHILIC GAUZE

- Size: 36x40 cm
- Threads: 12/8
- Minimum order: carton of 48 box of 12 pcs

Folded hydrophilic gauze 100% pure cotton, suitable for sterile medical dressings, to buffer and protect wounds.
It allows skin perspiration.
Individually packed and latex free.

Instructions and box in GB, FR, IT, ES, DE, GR.
Made in Italy.