

# ELASTICATED ADHESIVE BANDAGE 6 cm x 2.5 m non stretched



|               |                |
|---------------|----------------|
| Product Code  | 34783          |
| Unit of sale  | box of 10 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| RDM (NSIS) | 1848751       |
| CND        | M0304010301   |
| EAN/UPC    | 7630016124877 |
| PARAF      | 987474855     |

## Description

Adhesive bandage size 2.5 m x 6 cm  
(Length: 2.5 m non stretched, up to 4.5 m stretched)

Plaster bandages with hypoallergenic, synthetic adhesive containing zinc oxide, for permanent compression and functional dressings.  
Hypoallergenic, porous, air permeable with woven edges, can be sterilized by ethylene oxide.  
100% cotton, latex free.

MADE IN SWITZERLAND.