

# ELASTICATED ADHESIVE BANDAGE 8 cm x 2.5 m non stretched



Product Code	34784
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1848751
CND	M0304010301
EAN/UPC	7630016124846
PARAF	987474867

## Description

34784 Adhesive bandage size 2.5 m x 8 cm  
(Length: 2.5 m non stretched, up to 4.5 m stretched)

Plaster bandages with hypoallergenic, synthetic adhesive containing zinc oxide, for permanent compression and functional dressings. Hypoallergenic, porous, air permeable with woven edges, can be sterilized by ethylene oxide. 100% cotton, latex free.

MADE IN SWITZERLAND.