

ELASTICATED ADHESIVE BANDAGE 10 cm x 2.5 m non stretched



Product Code 34785

Unit of sale box of 10 pcs

Minimum order 1

Type Medical device

Class I

RDM (NSIS) 1848751

CND M0304010301

EAN/UPC 7630016124976

PARAF 987474842

Description

Adhesive bandage size 2.5 m x 10 cm

(Length: 2.5 m non stretched, up to 4.5 m stretched)

Plaster bandages with hypoallergenic, synthetic adhesive containing zinc oxide, for permanent compression and functional dressings. Hypoallergenic, porous, air permeable with woven edges, can be sterilized by ethylene oxide. 100% cotton, latex free.

MADE IN SWITZERLAND.