

DIGITAL SILICONE NURSE WATCH - blue



Product Code	20480
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279204803

Description

Digital nurse Watch

Colour: blue

Silicone nurse watch, provided with a safety pin.

Suitable to measure patient pulse, beat and sterilization procedure timing.

Conforming to procedures not allowing nurses to wear wrist watches.

- Watch case 35 mm
- Watch length 85 mm
- Manual: GB, FR, IT, ES