

PLASTER DISPENSER with 2 refill



Product Code	34826
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279348262

Description

PLASTER DISPENSER

Supplied with 2 refills (one each type):

- Code 34827 REFILL 42 TEXTILE PLASTERS 19 x 72 mm
- Code 34828 REFILL 36 TEXTILE PLASTERS 25 x 72 mm

Dispenser for adhesive plasters ready for use. It allows operator to extract plaster with one hand and apply directly on the wound.
Transparent cover protects plaster and shows residual quantities.
Lockable and with fixing wall kit included.

Technical Specifications

Size: 230 x 120 x 38 mm