

COTTON WOOL 500 g - Z folded



Product Code34835

Unit of salebox of 20 pcs

Minimum order1

TypeMedical device

ClassI

UK-REPNo

CH-REPNo

RDM (NSIS)2557394

CNDM010101

EAN/UPC8023279348354

PARAF986008530

Description

100% COTTON WOOL

- Cotton wool 500 g
- Type: Z folded
- Closing laces: No

Hydrophilic 100% cotton for sanitary and cosmetic use, hygiene, exudate absorption, non-injured skin disinfection.

Non-invasive Class I medical device.

Non sterile.

Multilingual pack: GB, FR, IT, DE, ES, FI, GR, HR, PL, PT, RO, SE.