



Product Code	34890
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2318584
CND	M04040502
GTIN	04054596954083

## Description

3M™ NU-GEL™

3M™ code: MNG415

3M™ Nu-Gel™ Hydrogel with Alginate

Minimum order: Box of 10

3M™ Nu-gel™ is a transparent, amorphous hydrogel containing sodium alginate to gently and effectively debride necrotic tissue and fibrinous slough.

It creates a moist wound healing environment which assists with natural autolytic debridement by donating moisture to the wound bed whilst the alginate component enhances its absorptive capabilities.

The gel can also be used to soften and hydrate eschar by facilitating rehydration of the wound.

Longer wear time.