

TEGADERM 3M - 6x7 cm for IV canulae - sterile



Product Code	34894
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2265722
CND	M04040601
EAN/UPC	4001895928734
PARAF	977814817

Description

3M™ TEGADERM™ FILM

- 3M Code: 1623W
- Product: Tegaderm™ - 6x7 cm (For I.V. cannule) - sterile
- Minimum Order: Box of 100

Transparent, sterile waterproof dressing.
Highly breathable.
Frame style for easy application also available in non-sterile roll.