

TEGADERM 3M + PAD - 5x7 cm - sterile



Product Code	34895
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2393770
CND	M04010102
EAN/UPC	8711428086881
PARAF	977814829

Description

3M medication products - Sterile transparent film dressing - 3582
 TEGADERM™ - PAD 5 x 7 cm - sterile
 Minimum order: 1 box of 50 pcs

Transparent, sterile waterproof dressing. Highly breathable.
 Frame style for easy application.
 With a non-adherent absorbent PAD.
 Ideal for surgical wounds.