

TEGADERM 3M - 10 cm x 10 m - non sterile



Product Code	34899
Unit of sale	1 roll
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2753742
CND	M0599
EAN/UPC	8023279348996
PARAF	977814868

Description

3M™ TEGADERM™ FILM

- 3M Code: 16004
- Product: Tegaderm™ - 10 cmx10 m non-sterile
- Minimum Order: 1 roll

Transparent, sterile waterproof dressing.
Highly breathable.
Frame style for easy application also available in non-sterile roll.