

FINGER GONIOMETER 6" - stainless steel



Product Code	24900
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279249002
---------	---------------

Description

FINGER GONIOMETER 6" - 15 cm

Finger goniometer measures range of motion of finger joints.
Linear scale in inches and centimeters on arms.
180° scales marked in 5° increments.

Multilanguage manual: GB, FR, IT, ES.

Approximate measurement of values.