

LUXIFLEX PLUS II LED LIGHT 50,000 Lux - trolley



Product Code	49021
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

Description

LUXIFLEX II PLUS on trolley - LED

Powerful LED light source that offers exceptional light quality. Flexible, height-adjustable and highly precise arm allows adjustment to any position. Compact head design and switch located at the base of the lamp. Class II.

Led Plus II with CRI 97%

Plus II has the high performance of the color Ra (97%) that depicts the exceptional quality of white light for an excellent distinction of fabrics and reproduction of all visible colours, naturally and precisely.

Manual in: GB, FR, IT, DE, ES, PT

Made in Spain.

Technical Specifications

Consumption (Watt) 12
Life of the bulbs (hours) 50,000
Light intensity at 50 cm (lux) 50,000
Light intensity at 1 m (lux) 17,000
Colour temperature $\pm 5\%$ at 1 m (K) 4,500°
Colour rendering index (Ra)% 97
Light field size at 50 cm (Ø mm) 150
External Ø of the reflector (mm) 65
Colour white
Norms EN 60601-1-2
Marks EU UE 2017/745
Weight (kg): -
Operating voltage V 100/240
Operating Hz 50/60