



Product Code	49035
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2570223
CND	Z129004
EAN/UPC	8023279490350

## Description

LUMINA LED LIGHT - trolley

Classic examination LED light with flexible arm.

Manual: GB, FR, IT, ES, PT, DE, PL, HU, SK, RO, BG, HR, SI, GR.

Flexible arm, extension 76 cm adjustable.  
Integrated handle on lamphead.  
360° degree free position with PVC coated, flexible metal arm.  
Brightness control.  
ABS base with brake.

## Technical Specifications

Consumption: 3 W  
Light intensity at 50/30 cm: 40,000/100,000 lux  
Colour temperature: 6,000 °K  
External Ø of reflector: 75 mm  
Colour rendering index (CRI): 90 Ra  
Light field Ø: 4 cm at 50 cm  
Operating voltage: 110-265 V, 50/60 Hz