

180° EXTRA LONG ARM GONIOMETER 14" - stainless steel



Product Code	24904
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279249040
---------	---------------

Description

180° EXTRA LONG ARMS GONIOMETER 14" - 35 cm

Extra long arms, useful for measuring all joints.
Arm has a linear scale marked in 1° increments.

Multilanguage manual: GB, FR, IT, ES.

Approximate measurement of values.