

HYRIDIA 7 LEDS LIGHT with metal spring arm - trolley



Product Code	49041
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2570223
CND	Z129004
EAN/UPC	8023279490411

Description

HYRIDIA 7 LEDS LIGHT - metal spring arm - trolley

7 leds 10W medical procedure lights on trolley with flexible soft arm or metal spring arm.

Digital brightness control.

Manual: GB, FR, IT, ES, PT, DE, PL, HU, SK, RO, BG, HR, SI, GR.

Handle on the lamphead.
Spring metal arm.
Detachable large base with brake.

Technical Specifications

- Consumption: 15 W
- Colour temperature: 4,800-5,500 °K
- External Ø of reflector: 89 mm
- Colour rendering index (CRI): 90 Ra
- Operating voltage: 110-265 V, 50/60 Hz
- Light intensity and light field Ø:
illumination 180,000 lux / Ø 9 cm / distance 30 cm
illumination 100,000 lux / Ø 12 cm / distance 50 cm
illumination 30,000 lux / Ø 18 cm / distance 100 cm