

LUMEX LED LIGHT with handle - trolley



Product Code	49045
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2570226
CND	Z129004
EAN/UPC	8023279490459

Description

LUMEX LED LIGHT with handle - trolley

Compact LED light featuring a modern design combined with innovative lighting technology and excellent manoeuvrability. The light is significantly more pleasant to the eyes and colours are rendered accurately. Ideal for examinations and minor surgery procedures at doctor's office, clinics and hospitals. Ergonomic and hygienic design with removable and sterilizable handle which allows movements.

Control system to manage 5 levels of brightness

Directional handle

Multilingual manual: GB, FR, IT, ES, PT, DE, PL, NL, CZ, HU, SE, DK, HR, RO, GR.

Technical Specifications

Light Intensity at 50/100 cm: 45,000 / 18,000 lux
Control System: 5 levels of light-touch brightness control and on/off
Lamphead: Ø 120 mm, 28 mm thickness
Diameter of illuminated field: 180-300 mm at 50-100 cm distance
Light source: 6 LEDs
LED life: 50,000 hours
Colour temperature: 3500-4500 K
Colour rendering index (CRI): 95
Power supply: 100-240V, 50/60 MHz
Power consumption: 6W
Protection: Class II