



Product Code	49047
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	8023279490473

Description

LED SPOT-ADJUSTED LIGHT - trolley

LED Spot examination lamp with multi-optical lens and light spot adjustable by inner mini light grating.

Additional features:

- uniform and pure light spot
- front optical lens
- inner light grating
- rear dual optical lens
- 100 cm super long 360° free position and PVC coated flexible metal arm

Manual: GB, FR, IT, ES, PT, DE, PL, GR.

Stainless steel handle.

Super long flexible soft arm, length 100 cm.

Light spot adjustable by circle on lamp head.

Digital brightness control.

ABS module base.

Technical Specifications

Bulbs/power: 1 LED / 5 W

Ø of lamp head: 25 mm

Illumination:

25,000 lux at 30 cm

10,000 lux at 50 cm

Light field Ø: max 70 mm, min 15 mm at 30 cm

Colour temperature: 4,800°-6,000° K

Lamphead: Ø 25 mm, length 125 mm

Power supply: 100-260 V, 50/60 Hz

Power cord: length 272 cm