



Product Code	49050
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	1189466
CND	Z12010701
EAN/UPC	4045396168251

Description

RI-MAGIC® HP LED - desk

High performance medical LED light for examination and procedures.

- high performance LED, 6,000 °K
- lifetime of the lamp up to 50,000 hours
- significantly reduced heating and energy loss
- focusing device at the front of the light head
- highly flexible, long special arm can be bent in all directions
- rheostat handle for infinitely variable regulation of light intensity.

Manual in: GB, FR, IT, ES, DE, RU.

Made in Germany.

Technical Specifications

- Distance: 1000 mm - 800 mm - 600 mm - 400 mm - 200 mm
- Lux max: 15,000 - 20,000 - 30,000 - 45,000 - 55,000
- Circle of light max Ø: 340 mm - 270 mm - 215 mm - 155 mm - 90 mm
- Circle of light min Ø: 40 mm - 35 mm - 30 mm - 25 mm - 20 mm
- Operating voltage: 100-240V, 50-60 Hz